

(E/Z)-Zotiraciclib citrate

Chemical Properties

CAS No. :

Formula: C₂₉H₃₂N₄O₈

Molecular Weight: 564.59

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	(E/Z)-Zotiraciclib ((E/Z)-TG02) citrate is a potent inhibitor of CDK2, JAK2 and FLT3 and can be used in cancer research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7712 mL	8.856 mL	17.712 mL
5 mM	0.3542 mL	1.7712 mL	3.5424 mL
10 mM	0.1771 mL	0.8856 mL	1.7712 mL
50 mM	0.0354 mL	0.1771 mL	0.3542 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481