

SPOP-IN-6b dihydrochloride

Chemical Properties

CAS No. :

Formula: C₂₈H₃₄Cl₂N₆O₃

Molecular Weight: 573.51

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	SPOP-IN-6b dihydrochloride is a potent inhibitor of spotted POZ protein (SPOP) (IC ₅₀ : 3.58 μM).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7436 mL	8.7182 mL	17.4365 mL
5 mM	0.3487 mL	1.7436 mL	3.4873 mL
10 mM	0.1744 mL	0.8718 mL	1.7436 mL
50 mM	0.0349 mL	0.1744 mL	0.3487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481