Data Sheet (Cat.No.T64103)



RET-IN-12

Chemical Properties


CAS No.: 2684252-55-7

Formula: C30H30F3N5O4

Molecular Weight: 581.59

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description RET-IN-12 is a RET inhibitor capable of acting on RET(WT) (IC50: 0.3 nM) and RET(V804M)

(IC50: 1 nM).

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7194 mL	8.5971 mL	17.1942 mL
5 mM	0.3439 mL	1.7194 mL	3.4388 mL
10 mM	0.1719 mL	0.8597 m ⊡	1.7194 mL
50 mM	0.0344 mL	0.1719 mL	0.3439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com