

Heparin Lithium salt

Chemical Properties

CAS No. : 9045-22-1
Formula: C₁₄H₂₅NO₂OS₃
Molecular Weight:
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Heparin Lithium salt

Biological Description

Description	Heparin Lithium salt is an anticoagulant that reversibly binds antithrombin III (ATIII). Heparin Lithium salt markedly blocks exosome-cell interactions.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481