

Laropiprant sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Laropiprant sodium is a potent and selective DP receptor antagonist with Ki values of 0.57 nM for TP receptors and 2.95 nM for DP receptors.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481