

Pacritinib hydrochloride

Chemical Properties


CAS No. : 1228923-43-0

Formula: C₂₈H₃₂N₄O₃.xClH

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Pacritinib hydrochloride is a potent inhibitor of wild-type JAK2 (IC ₅₀ : 23 nM) and the JAK2V617F mutant (IC ₅₀ : 19 nM). Pacritinib hydrochloride can be used in the study of acute myeloid leukaemia (AML) and myelofibrosis (MF).
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481