


2,3,5-Trimethylphenol

Chemical Properties

CAS No. :	697-82-5
Formula:	C ₉ H ₁₂ O
Molecular Weight:	136.19
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	2,3,5-Trimethylphenol (Isopseudocumenol) is a antiproliferative agent, pGI50= 0.36 to T. pyriformis.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 50 mg/mL (367.13 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.3427 mL	36.7134 mL	73.4268 mL
5 mM	1.4685 mL	7.3427 mL	14.6854 mL
10 mM	0.7343 mL	3.6713 mL	7.3427 mL
50 mM	0.1469 mL	0.7343 mL	1.4685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Schultz TW. Structure-toxicity relationships for benzenes evaluated with *Tetrahymena pyriformis*. Chem Res Toxicol. 1999 Dec;12(12):1262-7.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481