


3-pyridine toxoflavin

Chemical Properties

CAS No. :	32502-20-8
Formula:	C ₁₂ H ₁₀ N ₆ O ₂
Molecular Weight:	270.25
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	3-pyridine toxoflavin is a potent KDM4C inhibitor, IC ₅₀ =19nM.
Targets(IC ₅₀)	Histone Demethylase

Solubility Information

Solubility	DMSO: 45 mg/mL (166.51 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7003 mL	18.5014 mL	37.0028 mL
5 mM	0.7401 mL	3.7003 mL	7.4006 mL
10 mM	0.370 mL	1.8501 mL	3.7003 mL
50 mM	0.074 mL	0.370 mL	0.7401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Letfus, Vatroslav; Bioorganic & Medicinal Chemistry (2020), 28(1), 115128, CAPLUS and MEDLINE

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481