


Compound V007-6018

Chemical Properties

CAS No. : 859093-06-4
Formula: C₃₀H₂₄F₆N₄O₂S
Molecular Weight: 618.59
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Compound V007-6018 is a pharmaceutical compound used for drug screening.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 6.19 mg/mL (10 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6166 mL	8.0829 mL	16.1658 mL
5 mM	0.3233 mL	1.6166 mL	3.2332 mL
10 mM	0.1617 mL	0.8083 mL	1.6166 mL
50 mM	0.0323 mL	0.1617 mL	0.3233 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481