


Mifentidine

Chemical Properties

CAS No. :	83184-43-4
Formula:	C ₁₃ H ₁₆ N ₄
Molecular Weight:	228.29
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Mifentidine (DA 4577) is an orally available H ₂ receptor antagonist for the study of gastric emptying and experimental gastric and duodenal ulcers.
Targets(IC50)	Histamine Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3804 mL	21.902 mL	43.8039 mL
5 mM	0.8761 mL	4.3804 mL	8.7608 mL
10 mM	0.438 mL	2.1902 mL	4.3804 mL
50 mM	0.0876 mL	0.438 mL	0.8761 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Francalanza G, Franco C, Giachetti A. Mifentidine: evaluation of antiandrogen effects and kidney function studies. *Arzneimittelforschung*. 1985;35(1A):456-8.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481