


Clopentixol

Chemical Properties

CAS No. : 982-24-1
Formula: C₂₂H₂₅ClN₂O₂
Molecular Weight: 400.96
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Clopentixol is a thioxanthene with therapeutic actions similar to the phenothiazine antipsychotics. It is an antagonist at D1 and D2 dopamine receptors.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.494 mL	12.4701 mL	24.9401 mL
5 mM	0.4988 mL	2.494 mL	4.988 mL
10 mM	0.2494 mL	1.247 mL	2.494 mL
50 mM	0.0499 mL	0.2494 mL	0.4988 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481