


MI-63

Chemical Properties

CAS No. : 908027-21-4
Formula: C₂₉H₃₅Cl₂FN₄O₃
Molecular Weight: 577.52
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	MI-63 is a potent, specific, non-peptide inhibitor of the MDM2-p53 interaction.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7315 mL	8.6577 mL	17.3154 mL
5 mM	0.3463 mL	1.7315 mL	3.4631 mL
10 mM	0.1732 mL	0.8658 mL	1.7315 mL
50 mM	0.0346 mL	0.1732 mL	0.3463 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481