


AG-28262 besylate

Chemical Properties

CAS No. : 881688-70-6
Formula: C₂₈H₂₄N₄O₅S₃
Molecular Weight: 592.7
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AG-28262 besylate is a VEGFR-2 Inhibitor which may affect alanine aminotransferase gene expression and enzymatic activity in the liver.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6872 mL	8.436 mL	16.8719 mL
5 mM	0.3374 mL	1.6872 mL	3.3744 mL
10 mM	0.1687 mL	0.8436 mL	1.6872 mL
50 mM	0.0337 mL	0.1687 mL	0.3374 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481