## Data Sheet (Cat.No.T68546)



## PF-00337210

## **Chemical Properties**

CAS No.: 854514-88-8

Formula: C26H27N3O5

Molecular Weight: 461.51

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

## **Biological Description**

Description PF-00337210 is an orally available ATP-competitive inhibitor of the vascular endothelial

growth factor receptor type 2 (VEGFR2), with potential anti-angiogenesis and antineoplastic activities. Upon administration, the VEGFR2 tyrosine kinase inhibitor PF-00337210 selectively binds to VEGFR2 and prevents its phosphorylation which may result in an inhibition of migration, proliferation and survival of endothelial cells, microvessel formation, the inhibition of tumor cell proliferation, and may eventually cause tumor cell death. VEGFR2, a receptor tyrosine kinase, is frequently overexpressed by a variety of

tumor types.

**Preparing Stock Solutions** 

	<u> </u>			
	1mg	5mg	10mg	
1 mM	2.1668 mL	10.834 mL	21.668 mL	
5 mM	0.4334 mL	2.1668 mL	4.3336 mL	
10 mM	0.2167 mL	1.0834 mL	2.1668 mL	
50 mM	0.0433 mL	0.2167 mL	0.4334 mL	

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com