

AGN-204396

Chemical Properties


CAS No. : 847665-57-0

Formula: C32H44FN3O4

Molecular Weight: 553.71

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AGN-204396 is an antagonist that selectively blocks the activity of prostamides (prostaglandin-ethanolamides).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.806 mL	9.030 mL	18.060 mL
5 mM	0.3612 mL	1.806 mL	3.612 mL
10 mM	0.1806 mL	0.903 mL	1.806 mL
50 mM	0.0361 mL	0.1806 mL	0.3612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481