


Deterenol Free Base

Chemical Properties

CAS No. : 7376-66-1
Formula: C₁₁H₁₇NO₂
Molecular Weight: 195.26
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Deterenol Free Base is a beta-adrenoceptor agonist. It is an effective nonmydriatic and nonmiotic hypotensive agent, which can be used in antiglaucoma treatment.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1214 mL	25.6069 mL	51.2138 mL
5 mM	1.0243 mL	5.1214 mL	10.2428 mL
10 mM	0.5121 mL	2.5607 mL	5.1214 mL
50 mM	0.1024 mL	0.5121 mL	1.0243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481