


Wy 25093

Chemical Properties

CAS No. : 72807-01-3
Formula: C₂₂H₂₄N₄O₂
Molecular Weight: 376.45
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Wy 25093 is an inhibitor of the neuronal 5-hydroxytryptamine (5-HT) uptake process in vitro and in vivo, and may possess potential antidepressant activity.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6564 mL	13.282 mL	26.564 mL
5 mM	0.5313 mL	2.6564 mL	5.3128 mL
10 mM	0.2656 mL	1.3282 mL	2.6564 mL
50 mM	0.0531 mL	0.2656 mL	0.5313 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481