


LDT3 HCl

Chemical Properties

CAS No. : 693227-88-2
Formula: C₂₀H₂₈Cl₂N₂O₂
Molecular Weight: 399.35
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LDT3 HCl is a potent multi-target antagonists of both alpha1A- and alpha1D-adrenoceptors, and also of 5-HT1A receptors.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5041 mL	12.5203 mL	25.0407 mL
5 mM	0.5008 mL	2.5041 mL	5.0081 mL
10 mM	0.2504 mL	1.252 mL	2.5041 mL
50 mM	0.0501 mL	0.2504 mL	0.5008 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481