


Saringosterol

Chemical Properties

CAS No. : 6901-60-6
Formula: C₂₉H₄₈O₂
Molecular Weight: 428.69
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Saringosterol, also known as Saningosterol and Cystosphaerol, is a selective LXR β agonist. Saringosterol is racemic at the 24-position.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3327 mL	11.6634 mL	23.3269 mL
5 mM	0.4665 mL	2.3327 mL	4.6654 mL
10 mM	0.2333 mL	1.1663 mL	2.3327 mL
50 mM	0.0467 mL	0.2333 mL	0.4665 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481