


(-)-Adenophorine

Chemical Properties

CAS No. : 616215-43-1
Formula: C₈H₁₇NO₄
Molecular Weight: 191.22
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	(-)-Adenophorine is a moderate alpha-L-fucosidase inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2296 mL	26.1479 mL	52.2958 mL
5 mM	1.0459 mL	5.2296 mL	10.4592 mL
10 mM	0.523 mL	2.6148 mL	5.2296 mL
50 mM	0.1046 mL	0.523 mL	1.0459 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481