Data Sheet (Cat.No.T69000)



Prenalterol hydrochloride

Chemical Propert	ties
CAS No. :	61260-05-7
Formula:	C12H20ClNO3
Molecular Weight:	
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year

Biological Description	e e
Description	Prenalterol hydrochloride is a partial adrenergic agonist with functional beta 1-receptor specificity and inotropic effect. It is effective in the treatment of acute cardiac failure, postmyocardial infarction low-output syndrome, shock, and reducing orthostatic
	hypotension in the Shy-Drager syndrome.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8204 mL	19.1022 mL	38.2044 mL
5 mM 🛞	0.7641 mL	3.8204 mL	7.6409 mL
10 mM	0.382 mL	1.9102 mL	3.8204 mL
50 mM	0.0764 mL	0.382 mL	0.7641 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor • Natural Compounds • Compound Libraries • Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481