


Endosidin1

Chemical Properties

CAS No. : 57672-91-0
Formula: C₃₈H₅₀O₁₆
Molecular Weight: 762.802
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Endosidin1 is a natural vesicle trafficking inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.311 mL	6.5548 mL	13.1096 mL
5 mM	0.2622 mL	1.311 mL	2.6219 mL
10 mM	0.1311 mL	0.6555 mL	1.311 mL
50 mM	0.0262 mL	0.1311 mL	0.2622 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481