


Chondrocurine

Chemical Properties

CAS No. : 477-58-7
Formula: C₃₆H₃₈N₂O₆
Molecular Weight: 594.7
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Chondrocurine is a biochemical that inhibits quick death of mice treated with bacillus Calmetta Guerin & elicited with lipopolysaccharide.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6815 mL	8.4076 mL	16.8152 mL
5 mM	0.3363 mL	1.6815 mL	3.363 mL
10 mM	0.1682 mL	0.8408 mL	1.6815 mL
50 mM	0.0336 mL	0.1682 mL	0.3363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481