


P-1946

Chemical Properties

CAS No. : 449806-40-0
Formula: C₂₉H₄₁N₅O₆S
Molecular Weight: 587.73
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	P-1946 is an HIV-1 protease inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7015 mL	8.5073 mL	17.0146 mL
5 mM	0.3403 mL	1.7015 mL	3.4029 mL
10 mM	0.1701 mL	0.8507 mL	1.7015 mL
50 mM	0.034 mL	0.1701 mL	0.3403 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481