


Ascofuranone

Chemical Properties

CAS No. : 38462-04-3
Formula: C₂₃H₂₉ClO₅
Molecular Weight: 420.93
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Ascofuranone is an inhibitor of the ubiquinol oxidase activity of <i>Trypanosoma brucei</i> mitochondrial alternative oxidase (TAO), as well as an inhibitor of HsDHODH.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3757 mL	11.8785 mL	23.7569 mL
5 mM	0.4751 mL	2.3757 mL	4.7514 mL
10 mM	0.2376 mL	1.1878 mL	2.3757 mL
50 mM	0.0475 mL	0.2376 mL	0.4751 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481