


Pyrvinium Iodide

Chemical Properties

CAS No. : 35648-29-4
Formula: C₂₆H₂₈IN₃
Molecular Weight: 509.4355
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Pyrvinium Iodide is a casein kinase 1α-dependent inhibitor of Wnt/β-Catenin.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9629 mL	9.8147 mL	19.6294 mL
5 mM	0.3926 mL	1.9629 mL	3.9259 mL
10 mM	0.1963 mL	0.9815 mL	1.9629 mL
50 mM	0.0393 mL	0.1963 mL	0.3926 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481