


CNS-5788

Chemical Properties

CAS No. : 342047-49-8
Formula: C₁₆H₁₈ClN₃OS₂
Molecular Weight: 367.92
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	CNS-5788 is a potential ischemia-selective N-methyl-D-aspartate (NMDA) receptor ion-channel blocker
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.718 mL	13.5899 mL	27.1798 mL
5 mM	0.5436 mL	2.718 mL	5.436 mL
10 mM	0.2718 mL	1.359 mL	2.718 mL
50 mM	0.0544 mL	0.2718 mL	0.5436 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481