


CGC 11150

Chemical Properties

CAS No. : 304911-08-8
Formula: C₄₀H₁₀₀Cl₁₀N₁₀
Molecular Weight: 1075.81
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	CGC 11150 is a polyamine analogue inhibit tumor growth in vitro and in vivo.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9295 mL	4.6477 mL	9.2953 mL
5 mM	0.1859 mL	0.9295 mL	1.8591 mL
10 mM	0.093 mL	0.4648 mL	0.9295 mL
50 mM	0.0186 mL	0.093 mL	0.1859 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481