

R 101933

Chemical Properties


CAS No. : 278798-78-0

Formula: C37H38N4O3

Molecular Weight: 586.72

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description

R 101933 is a P-glycoprotein (P-gp) inhibitor.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7044 mL	8.522 mL	17.0439 mL
5 mM	0.3409 mL	1.7044 mL	3.4088 mL
10 mM	0.1704 mL	0.8522 mL	1.7044 mL
50 mM	0.0341 mL	0.1704 mL	0.3409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:36 Washington Street,Wellesley Hills,MA 02481