


8RK59

Chemical Properties

CAS No. : 2705841-53-6
Formula: C₃₃H₃₉BF₂N₁₀O₂S
Molecular Weight: 688.6
Appearance: no data available
Storage: keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	8RK59 is a PARK7 inhibitor.
-------------	-----------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4522 mL	7.2611 mL	14.5222 mL
5 mM	0.2904 mL	1.4522 mL	2.9044 mL
10 mM	0.1452 mL	0.7261 mL	1.4522 mL
50 mM	0.029 mL	0.1452 mL	0.2904 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481