


VUF16839

Chemical Properties

CAS No. : 2387674-52-2
Formula: C₁₀H₁₇N₅
Molecular Weight: 207.28
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	VUF16839 is a high-affinity non-imidazole histamine H3 receptor agonist.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8244 mL	24.122 mL	48.2439 mL
5 mM	0.9649 mL	4.8244 mL	9.6488 mL
10 mM	0.4824 mL	2.4122 mL	4.8244 mL
50 mM	0.0965 mL	0.4824 mL	0.9649 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481