


BI-0314

Chemical Properties

CAS No. : 2244560-46-9
Formula: C₁₃H₁₇F₃N₄
Molecular Weight: 286.3
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	BI-0314 is a selective allosteric STEP activator.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4928 mL	17.4642 mL	34.9284 mL
5 mM	0.6986 mL	3.4928 mL	6.9857 mL
10 mM	0.3493 mL	1.7464 mL	3.4928 mL
50 mM	0.0699 mL	0.3493 mL	0.6986 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481