


HMR-1566

Chemical Properties

CAS No. : 223749-45-9
Formula: C₁₈H₂₆F₃NO₅S
Molecular Weight: 425.46
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	HMR-1566 is a potent and selective Iks channel blocker.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3504 mL	11.752 mL	23.504 mL
5 mM	0.4701 mL	2.3504 mL	4.7008 mL
10 mM	0.235 mL	1.1752 mL	2.3504 mL
50 mM	0.047 mL	0.235 mL	0.4701 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481