


YW1159

Chemical Properties

CAS No. : 2131223-69-1
Formula: C₁₉H₁₄FN₅O
Molecular Weight: 347.35
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	YW1159 is an inhibitor of Wnt signaling.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8789 mL	14.3947 mL	28.7894 mL
5 mM	0.5758 mL	2.8789 mL	5.7579 mL
10 mM	0.2879 mL	1.4395 mL	2.8789 mL
50 mM	0.0576 mL	0.2879 mL	0.5758 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481