


TargetMol®
YOUR TARGET MOLECULES

Chemical Properties

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Description	PF-06862309 is a PI3K inhibitor
-------------	---------------------------------

	1mg	5mg	10mg
1 mM	1.7288 mL	8.6439 mL	17.2879 mL
5 mM	0.3458 mL	1.7288 mL	3.4576 mL
10 mM	0.1729 mL	0.8644 mL	1.7288 mL
50 mM	0.0346 mL	0.1729 mL	0.3458 mL

www.targetmol.com