


Mesaconitine

Chemical Properties

CAS No. :	2752-64-9
Formula:	C33H45NO11
Molecular Weight:	631.71
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	1. Mesaconitine has antiinflammatory activity. 2. Mesaconitine has antinociceptive activity through cyclic AMP and stimulation of the central β -adrenergic system.
Targets(IC50)	TNF

Solubility Information

Solubility	DMSO: 11 mg/mL (17.41 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.583 mL	7.915 mL	15.830 mL
5 mM	0.3166 mL	1.583 mL	3.166 mL
10 mM	0.1583 mL	0.7915 mL	1.583 mL
50 mM	0.0317 mL	0.1583 mL	0.3166 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ameri A . Inhibition of stimulus-triggered and spontaneous epileptiform activity in rat hippocampal slices by the Aconitum alkaloid mesaconitine.[J]. European Journal of Pharmacology, 1998, 342(2-3):183.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481