


Ginsenoside Rk3

Chemical Properties

CAS No. :	364779-15-7
Formula:	C ₃₆ H ₆₀ O ₈
Molecular Weight:	620.86
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	1. Ginsenoside Rk3 and Rh4 could have a role in treating inflammatory diseases. 2. Ginsenoside Rk3 is often used as a major ingredient of the compound preparation for ischemic heart diseases.
Targets(IC50)	Others,NF-κB

Solubility Information

Solubility	DMSO: 50 mg/mL (80.53 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6107 mL	8.0533 mL	16.1067 mL
5 mM	0.3221 mL	1.6107 mL	3.2213 mL
10 mM	0.1611 mL	0.8053 mL	1.6107 mL
50 mM	0.0322 mL	0.1611 mL	0.3221 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Sun J , Sun G , Meng X , et al. Ginsenoside RK3 Prevents Hypoxia-Reoxygenation Induced Apoptosis in H9c2 Cardiomyocytes via AKT and MAPK Pathway[J]. Evidence-Based Complementary and Alternative Medicine, 2013,

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481