Data Sheet (Cat.No.T70015)



GSK563

Chemical Properties


CAS No.: 1983928-04-6

Formula: C16H22O2

Molecular Weight: 246.34

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	GSK563 is a potent GVIA iPLA2 inhibitor (XI(50) 0.0000021, IC50 1 nM). GSK563 is 22 000
	times more active against GVIA iPLA2 than GIVA cPLA2. It was found to reduce β-cell
	apoptosis induced by proinflammatory cytokines, raising the possibility that it can be
, and the second	beneficial in countering autoimmune diseases, such as type 1 diabetes.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0594 mL	20.2972 mL	40.5943 mL
5 mM	0.8119 mL	4.0594 mL	8.1189 mL
10 mM	0.4059 mL	2.0297 mL	4.0594 mL
50 mM	0.0812 mL	0.4059 mL	0.8119 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com