


ER-38925

Chemical Properties

CAS No. : 196517-43-8
Formula: C₂₁H₁₇NO₃
Molecular Weight: 331.36
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	ER-38925 is a retinoic acid receptor (RAR) subtype α -selective agonist
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0179 mL	15.0893 mL	30.1787 mL
5 mM	0.6036 mL	3.0179 mL	6.0357 mL
10 mM	0.3018 mL	1.5089 mL	3.0179 mL
50 mM	0.0604 mL	0.3018 mL	0.6036 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481