## Data Sheet (Cat.No.T70065)



## JNJ-61803534

Chemical Proper	ties
CAS No. :	1917306-14-9
Formula:	C23H23Cl2F6N3O4S
Molecular Weight:	
Appearance:	no data available
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year

<b>Biological Description</b>	e	
Description	JNJ-61803534 is a potent, orally active RORyt inverse agonist, demonstrating anti- inflammatory activity through inhibition of IL-17A production in human CD4+ T cells under Th17 differentiation conditions, with an IC50 of 9.6 nM.	

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6067 mL	8.0333 mL	16.0666 mL
5 mM	0.3213 mL	1.6067 mL	3.2133 mL
10 mM	0.1607 mL	0.8033 mL	1.6067 mL
50 mM	0.0321 mL	0.1607 mL	0.3213 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481