


XM 412

Chemical Properties

CAS No. : 177932-89-7
Formula: C₃₅H₄₄N₄O₉S₂
Molecular Weight: 728.876
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	XM 412, also known as DMP 450, is an HIV protease inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.372 mL	6.8598 mL	13.7197 mL
5 mM	0.2744 mL	1.372 mL	2.7439 mL
10 mM	0.1372 mL	0.686 mL	1.372 mL
50 mM	0.0274 mL	0.1372 mL	0.2744 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481