


Epolactaene

Chemical Properties

CAS No. :	167782-17-4
Formula:	C ₂₁ H ₂₇ N ₀ O ₆
Molecular Weight:	389.44
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Epolactaene is a neuritogenic compound in human neuroblastoma cells which selectively inhibit the activities of mammalian DNA polymerases and human DNA topoisomerase II.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5678 mL	12.8389 mL	25.6779 mL
5 mM	0.5136 mL	2.5678 mL	5.1356 mL
10 mM	0.2568 mL	1.2839 mL	2.5678 mL
50 mM	0.0514 mL	0.2568 mL	0.5136 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481