


FTI-249

Chemical Properties

CAS No. : 161721-65-9
Formula: C₁₅H₂₃N₃O₃S₂
Molecular Weight: 357.49
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	FTI-249 is a Farnesyltransferase inhibitor, which potently inhibited FTase (IC ₅₀ = 100Å–200 nM).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7973 mL	13.9864 mL	27.9728 mL
5 mM	0.5595 mL	2.7973 mL	5.5946 mL
10 mM	0.2797 mL	1.3986 mL	2.7973 mL
50 mM	0.0559 mL	0.2797 mL	0.5595 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481