


QQN-00358

Chemical Properties

CAS No. : 1498300-35-8
Formula: C₂₆H₂₄F₃N₃O₄
Molecular Weight: 499.48
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	QQN-00358 is a novel potent and selective tankyrase (TNKS) inhibitor.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0021 mL	10.0104 mL	20.0208 mL
5 mM	0.4004 mL	2.0021 mL	4.0042 mL
10 mM	0.2002 mL	1.001 mL	2.0021 mL
50 mM	0.040 mL	0.2002 mL	0.4004 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481