


NU6310

Chemical Properties

CAS No. : 1453215-59-2
Formula: C₂₀H₂₅N₅O₃S
Molecular Weight: 415.51
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	NU6310 is a potent DCK2 inhibitor.
-------------	------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4067 mL	12.0334 mL	24.0668 mL
5 mM	0.4813 mL	2.4067 mL	4.8134 mL
10 mM	0.2407 mL	1.2033 mL	2.4067 mL
50 mM	0.0481 mL	0.2407 mL	0.4813 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481