


RB 104

Chemical Properties

CAS No. : 145253-48-1
Formula: C₂₁H₂₃IN₂O₆
Molecular Weight: 526.32
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	RB 104 is a potent inhibitor of EC 3.4.24.11 that is used to detect nanogram quantities of the enzyme by inhibitor gel electrophoresis.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.900 mL	9.4999 mL	18.9998 mL
5 mM	0.380 mL	1.900 mL	3.800 mL
10 mM	0.190 mL	0.950 mL	1.900 mL
50 mM	0.038 mL	0.190 mL	0.380 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481