

ALRN-6924

Chemical Properties


CAS No. : 1451199-80-6

Formula: C₉₃H₁₃₈N₂₀O₂₁

Molecular Weight: 1872.21

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	ALRN-6924 is a novel dual MDM4 and MDM2 inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5341 mL	2.6706 mL	5.3413 mL
5 mM	0.1068 mL	0.5341 mL	1.0683 mL
10 mM	0.0534 mL	0.2671 mL	0.5341 mL
50 mM	0.0107 mL	0.0534 mL	0.1068 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 36 Washington Street, Wellesley Hills, MA 02481