


S 12340

Chemical Properties

CAS No. : 144754-35-8
Formula: C₂₄H₃₈N₂O₃S
Molecular Weight: 434.64
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	S 12340 is a potent inhibitor of the oxidative modification of low-density lipoprotein.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3008 mL	11.5038 mL	23.0075 mL
5 mM	0.4602 mL	2.3008 mL	4.6015 mL
10 mM	0.2301 mL	1.1504 mL	2.3008 mL
50 mM	0.046 mL	0.2301 mL	0.4602 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481