

TargetMol®
YOUR TARGET MOLECULES

Chemical Properties

Formula: $C_{28}H_{44}O_3$

CCCCCCCCCCCCCCCCC(CCC)c1ccc(cc1)C(=O)/C=C/C=C

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Description	OBBA is a phospholipase A2 antagonist.
-------------	--

	1mg	5mg	10mg
1 mM	2.3329 mL	11.6645 mL	23.3291 mL
5 mM	0.4666 mL	2.3329 mL	4.6658 mL
10 mM	0.2333 mL	1.1665 mL	2.3329 mL
50 mM	0.0467 mL	0.2333 mL	0.4666 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481