


AGN-190121

Chemical Properties

CAS No. : 132032-67-8
Formula: C₂₀H₂₂O₂
Molecular Weight: 294.39
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AGN-190121 is a retinoic acid receptor (RAR) agonist. RAR and Retinoid X Receptor (RXR) ligands can act synergistically to induce hypertriglyceridemia through distinct mechanisms of action.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3969 mL	16.9843 mL	33.9685 mL
5 mM	0.6794 mL	3.3969 mL	6.7937 mL
10 mM	0.3397 mL	1.6984 mL	3.3969 mL
50 mM	0.0679 mL	0.3397 mL	0.6794 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481